

SWZ-1





SWP-2




Application: This series cabinet is the essential equipment to carry out biotechnology research and experimental, widely used in medical, college, farming and other industries Provide air purification for the working area.

Features:

Two models type, horizontal airflow and vertical airflow, According to using function it has desktop and standing type, single people single side (two sides) and two people single side (two sides), warranty 2 years.



Low energy consumption centrifugal fan and reasonable air duct design are adopted to make the air flow more uniform and stable. The operating table is mirror stainless steel for corrosion resistant cleaning.



High air pressure fan

The control system adopts microcomputer integration, adjustable wind speed, UV sterilization and lighting belt interlocking functions to effectively protect the operator;




Primary and high-efficiency filters are provided to make the filtering effect meet the user's purification level requirements;



- 1.The frame and cabinet adopt separate assembly design, convenient transportation, reduce transportation cost, and then when using Foot frame assembly.
- 2.Ergonomic tilt design for more comfortable operation, with folding glass front door
When not in use, it can be fully covered to reduce dust into the studio.

Model	Wind speed(m/s)	Illumination (LX)	Noise	Voltage (V)	Power(KW)	Inner size(mm) (H×W×D)	External size(mm) (H×W×D)	Packing size(mm) (H×W×D)	Net weight(KG)	Shipping weight(KG)
SW-650 (desktop)	0.25~0.45	≥300	≤62dB	220	0.15	480×650×55 0	680×660×5 80	700×720×6 40	40	43
SW-800 (desktop)					0.15	480×800×55 0	680×810×5 80	700×880×6 40	46	50
SWP-1(single person horizontal)					0.3	525×800×5 30	1515×860× 680	1200×960× 760	80	85
SWP-2 (double person horizontal)					0.5	525×1300× 530	1515×1360 ×680	1200×1480 ×760	120	127
SWZ-1(single person vertical)					0.3	510×800×6 50	1515×860× 680	1200×960× 760	80	85
SWZ-2 (double person vertical)					0.5	510×1300× 650	1515×1380 ×680	1200×1480 ×760	120	127

SWZ-1



SWZ-2



UV germicidal lamp



Lighting tube



control panel



Internal structure



Airflow screen



Stainless steel table



Move the wheel

Comparison	JHP horizontal	JHZ vertical	SW (Benchtop)	SWP horizontal	SWZ vertical
	Professional version		Standard version		
Microcomputer control					
LCD screen	●Filter failure reminder	●Filter failure reminder	●	●	●
Brand accessories					
ST bi type	●	●			
Philip chip	●	●			
Taiwan optical coupler	●	●			
Worktable					
SUS304 SS	●	●			
201 SS			●	●	●
Ergonomic design					
Front tilting design	●	●	●	●	●
Turned out design	●	●			
Positioned front door	●	●			
Folding front door			●	●	●
Centrifugal fan					
High temperature resistant bearing	●	●			
Hepa filter	●	●	●	●	●
Laminar flow					
Horizontal airflow	●			●	
Vertical airflow		●	●		●
Water and gas outfit, electric outfit	●	●			
Fan speed adjustable	●	●	●	●	●
CLEANFLOW technology	●	●			
Cleanliness rate	ISO Class5				
Interlocking lamps	●	●	●	●	●
Detached support	●	●		●	●
warranty	3 years		2 years		

**Manufacture of drying oven, lab incubator, climate chamber,
laminar flow cabinet, biological safety cabinet(OEM,ODM)**

KENTON APPARATUS LIMITED

KENTON Apparatus Limited was found in 1999, production plant was set in Guangzhou, China. It is mainly engaged in production, sales and export of laboratory thermostatic equipment and biological scientific equipment, and has provided high-quality products and services for professional laboratories and scientific research institutions at domestic and abroad for many years. We rely on an international professional R&D team to keep the design and technology of products in step with the world and constantly seek breakthroughs. Since its establishment, the company has been focusing on thermostatic equipment products in related fields. Now it is a high-tech enterprise integrating science, industry and trade. Its business covers China, Southeast Asia, America, Europe, Australia and other countries. At the same time, based on the business philosophy of "walking with scientific researchers", and taking advantage of the company's geographical advantages, the company has carried out extensive communication and exchange with scientific researchers in different regions for many years, and carried out long-term follow-up services to users, constantly improving the function and quality of products, so that products can better meet the requirements of users.

<https://www.kentonchina.com>

E-mail: kenton@kentonchina.com

**Adds: No.2, Hualong Rd, Donghua industrial zone,
Renhe town, Baiyun district, Guangzhou China**

Tel: 86-020-36246586, 36247961, 36246650, 36246651

